Seplife® 6AG Phe/HL/90



1. Description

Seplife® 6AG Phe/HL/90 is aromatic hydrophobic agarose resin designed for hydrophobic interaction chromatography (HIC) resin for capture and intermediate purification of biomolecules.

- The aromatic HIC resins is based on Seplife® 6AG/90 base matrix covalently functionalized with chemically-stable ether linkages.
- The agarose backbone ensures minimal levels of non-specific binding.
- The binding of biomolecules on Seplife® 6AG Phe/HL/90 is a HIC resin based on 4% cross-linked agarose and has a large particle size (45-165 micron) is performed at high salt concentration followed by elution with a low salt concentration solution.
- High stability to CIP (cleaning in place) up to 1M NaOH.
- Regulatory Support File (RSF) is available for Seplife® 6AG Phe/HL/90.

Seplife® 6AG Phe/HL/90 is a HIC resin based on highly cross-linked agarose (6%) and has a particle size range 45-165 micron.

2. Properties

Product	Seplife® 6AG Phe/HL/90	
Appearance	White spherical beads	
Туре	Strong hydrophobic chromatography media	
Matrix	6% cross-linked agarose	
Ligand	Phenyl	
Particle size range (μm)	45-165	
pH stability	3-13 (operational), 2-14 (CIP)	
Chemical Stability	Stable to commonly used aqueous buffers, 1.0 M NaOH , 3 M ammonium sulphate , 70% ethanol, 30% isopropanol, 0.5% SDS, 6 M guanidine hydrochloride, 8 M urea, 10% ethylene glycol	
Flow rate* (cm/h)	Max 750 cm/h	
Shipped as	20% ethanol slurry	

^{*}Testing conditions: Chromatography column 16mm \times 400mm; column bed height 25cm; temperature 25° C; mobile phase water.



Seplife® 6AG Phe/HL/90



3. Instructions

3.1 Column packing

Column loading should be performed in accordance with standard operating procedures. It is important to ensure that each material is at its working temperature, and when possible, the media should be degassed before column packing.

3.2 Equilibration

Equilibrate the column with an equilibration buffer solution of 2 to 5 times the volume of the column bed until the conductivity and pH of the effluent are completely consistent with those of the sample feeding buffer.

3.3 Sample feeding

- 1. Samples are prepared with an equilibration solution, and cloudy samples need to be centrifuged and filtered before loading.
- 2. Typically, the resin is used in bind-elute mode: let the target product bind to the column, wash off the impurities with the equilibration solution, and then choose an eluent to wash off the target product.
- 3. The extent to which the media adsorbs sample components depends on the hydrophobic nature of the sample, the ionic strength of the mobile phase, and the temperature. When the salt concentration is high, or the temperature is high, or the sample components are highly hydrophobic, the media will adsorb the components firmly.

3.4 Elution

The adsorbed biomolecules can be eluted by reducing the salt concentration. The elution can be enhanced by adding surfactants or organic solvents, the most commonly used is a low-salt concentration buffer, such as 0.02-0.05 mol/L PBS.

3.5 Regeneration

- 1. First wash with 3 to 5 CV of distilled water at the operating flow rate,. Then wash with 3 to 5 CV of the equilibration solution.
- 2. In case of inactivated proteins or lipids that cannot be washed away during regeneration, they can be removed by cleaning-in-place.

3.6 Cleaning-In-Place (CIP)

1. For proteins bound by ionic bonds, 0.5~1 BV of 2M NaCl can be used to remove them.



Seplife® 6AG Phe/HL/90



- 2. For precipitated proteins, hydrophobically bound proteins or lipids, first wash with 1 BV of 0.1M NaOH. Then wash with the equilibration buffer until the pH is neutral.
- 3. For proteins and lipids with strong hydrophobic binding, wash with 4-10 BV of 70% ethanol or 30% isopropanol.

It is important to note that the concentration of the organic solvent should gradually increase to avoid bubble formation.

4. Storage

Sealed and stored at 4~30°C (preservation solution is 20% ethanol) in a ventilated, dry and clean place, do not freeze.

5. Transportation

Avoid sunlight, rain, and heavy pressure during transportation. It is strictly forbidden to transport with toxic and hazardous materials.

6. Precaution

- **6.1** The order of the binding strength of typical ligands and proteins from strong to weak is the following: phenyl, octyl, butyl, isopropyl, ether. If the sample is tightly combined with the media, it can be replaced with a weaker media for easier elution.
- **6.2** In HIC, the hydrophobic chromatography media and hydrophobic ligands have a great influence on selectivity. Also, the influence of the following factors cannot be ignored: solubility of the sample, and the purification scale.
- **6.3** Since strong hydrophobic proteins and hydrophobic ligands are very tightly combined, when screening hydrophobic media, if the sample is a strong hydrophobic protein, start with low hydrophobic media and choose a media that obtains the highest resolution and capacity at low salt concentrations.
- **6.4** The sample and chromatography media Seplife®6AG Phe/HL/90must be thoroughly equilibrated with equilibration buffer before column chromatography can be performed.
- **6.5** Optimized separation can be obtained during gradient elution when using 1/5 of the total binding capacity of the chromatographic column.
- **6.6** The loaded column bed must have a flat surface, with no channel flow or air bubbles, otherwise it should be reloaded.



Seplife® 6AG Phe/HL/90



- **6.7** For large-scale purification and capture of proteins, step elution can reduce separation time and buffer consumption.
- **6.8** During the elution process, the flow rate should be strictly controlled.
- 6.9 During sample loading and the entire elution process, prevent the column surface from drying out.

7. Ordering information

Product Name	References	Pack Size
Seplife® 6AG Phe/HL/90	A3023202	25ml
	A3023203	100ml
	A3023204	500ml
	A3023205	1L
	A3023206	5L
	A3023207	10L

Production date: See label

Expiry date: 5 years, under proper storage conditions

Manufacturer: Sunresin New Materials Co. Ltd.

Add:No. 135, Jinye Rd, Xi'an Hi-tech Industrial Development Z one, Shaanxi, 710076, China

www.seplite.com www.sunresin.com
E-mail: info.lifescience@sunresin.com

All information set forth herein is for informational purposes only. This information is general descriptive(introductory) information of SUNRESIN and its related products, technologies and services. Neither shall constitute the guarantee of SUNRESIN and its affiliates to products, technologies and services in specific fields and specific application conditions results, unless otherwise expressly noted. SUNRESIN and its affiliates assumes no obligation or liability for the information in this document. Customer is responsible for judging whether the information is appropriate for Customer's concrete demand and are obliged to understand whether the use of these products, technologies and services is permitted by the laws and regulations of their countries and relevant regions. Unless expressly stated, no freedom from infringement of use any patent or trademark or intellectual property rights owned by SUNRESIN or its affiliated companies under this document is to be inferred.

