

## 1. Description

Seplife<sup>®</sup> LXPM SP5504M is an acrylic polymer for ion exchange chromatography characterized by higher hydrophilicity than styrene/DVB polymers for use in capture and polishing purification of medium size biomolecules.

- The surface of the polymer is modified to increase hydrophilicity followed by binding with strong acid ion exchange groups
- Narrow particle size range ideal for high resolution
- The polymer has a high porosity (1000Å) for biomolecules to bind without steric hindrance but high binding capacity even at high flow rates.
- High chemical stability to common aqueous solutions used in IEX chromatography, insoluble in most organic solvents.
- Regulatory Support File (RSF) is available

Seplife<sup>®</sup> LXPM SP5504M is a polymeric resin for ion exchange chromatography based on methacrylate functionalized with sulfopropyl (SP) with a uniform particle size (50-100 micron).

Product	Seplife® LXPM SP5504M	
Appearance	White to light yellow spherical beads	
Туре	Strong acid cation - Sulfopropyl	
Matrix	Polymethyl methacrylate	
lon exchange capacity (mmol/ml)	0.05-0.11	
pH ligand fully charged	Negatively Charged at pH>2.5	
Particle size range (micron)	50-100	
Typical pore size (Å)	1000	
pH stability	3-12 (operational), 3-14 (CIP)	
	All commonly used aqueous buffers: 1 M NaOH, 0.1 M acetic acid,	
Chemical Stability	8 M urea, 70% ethanol, 6 M guanidine hydrochloride, 2% benzyl	
	alcohol	
Flow rate* (cm/h)	500-600	
10% Dynamic binding capacity (mg /ml)**	≥40	
Shipped as	Slurry in 20% ethanol solution	

## 2.Properties

\*Testing conditions: Column size: 26x200mm, 0.3MPa, water

\*\*Report testing conditions: 20mM NaAc-HAc, pH 5.0, column size 8\*100, column height 10 cm; pH 5.0





### 3. Instructions

### 3.1 Column packing

The best column packing results are obtained with the resin slurry in 0.5M NaCl at a concentration of 60%-70%.

- Mix the resins to form a homogenate and measure the desired mass or volume, about 1.2 times the column volume.
- Replace 20% ethanol with 0.5 M NaCl solution and equilibrate overnight.
- Before loading the column, adjust the slurry concentration to 65-70% with 0.5 M NaCl solution; pour the resin slurry into the chromatography column in one go, allow to sediment by gravity and mark the bed height after sedimentation.
- Load the distributor and adjust the height so that the compression coefficient is 1.05-1.10; then pump 0.5M NaCl through the column and stabilize the bed with 1.5 to 2 times the highest working flow rate. (2 5 bed volumes (BV).
- Efficiency and asymmetry determinations can be performed as described below.

### 3.2 Column Efficiency Evaluation

The test method for column efficiency of the chromatography columns is as follows:

Mobile phase: 0.5 M NaCl solution

Linear flow rate: 100 cm/h

Sample: 2 M NaCl solution

Loading volume: 0.5-1% of column volume

Detection: Conductivity detector

Acceptable column asymmetry factor is in the range 0.8-1.8.

### 3.3 Rinsing

The loaded columns should be rinsed with at least 5 BV of deionized water.

#### 3.4 Equilibration

Equilibrate the column with an appropriate 5-10 column volume buffer until the conductivity and pH of the effluent remain unchanged. The specific buffer system should be screened and optimized according to the stability and isoelectric point of the target protein and the type of ion exchange medium.

#### 3.5 Sample feeding

The solid sample can be prepared by dissolving in the equilibrium solution; the low concentration sample solution can be dialyzed in the equilibrium buffer; the high concentration sample solution can be diluted with the equilibrium solution. To avoid clogging of the column, samples should be processed by centrifugation or membrane filtration. The feed amount is calculated according to the capacity of the resin and the content of the target protein in the feed solution. Before loading, ensure that the sample buffer is as consistent as possible with the equilibration solution.





### 3.6 Elution

After loading the sample, continue rinsing with equilibration buffer until the baseline is stable. Typically, the molecules retained on IEX chromatographic media are eluted by increasing the salt concentration or changing the pH to ensure repulsion of the charged species.

### 3.7 Regeneration

After each chromatography cycle, the media can be regenerated with 0.5-2 M NaCl to remove the strongly bound protein.

### 3.8 Cleaning-in-place (CIP)

In order to maintain the performance of the chromatographic column, if there are proteins or other impurities that cannot be effectively removed during the regeneration process, a CIP step should be performed. Up-flow CIP may increase the efficiency of the process. The recommended operation steps are as follows:

1) For precipitated, hydrophobically bound proteins or lipoproteins, wash with 0.2-0.5 M NaOH (contact time 1-2 h), followed by rinsing with equilibration solution (approx. 5 BV) and deionized water (approx. 3 BV);

2) For proteins retained by hydrophobic bonds, lipoproteins and lipids, wash with 50% ethanol or 30% isopropanol (approx. 5 BV, contact time 0.5-1 h), and rinse with deionized water (approx. 5 BV).

3) Other options: media can be cleaned with alkaline or acidic solutions containing non-ionic surfactants, such as 0.1-0.5% Triton X-100 in 0.1 M acetic acid for 1-2 h, and rinse with 50% ethanol (approx. 5 BV). To remove the detergent, rinse with water (approx. 5 BV). When using high-concentration organic solvents, it is recommended to gradually increase the solvent concentration in the mobile phase in order to avoid bubbling.

## 4. Storage

Chromatography resins that are not for immediate use should be stored in 20% ethanol at 4-30 °C. The loaded column should be stored in a buffer containing 20% ethanol (pH 7.0). Clean the ethanol after packing the column and before use.

### 5. Transportation

Avoid sunlight, rain, and heavy pressure during transportation. It is strictly forbidden to transport with toxic and hazardous materials.





## 6. Ordering information

Product Name	References	Pack Size
Seplife® LXPM SP5504M	PM50124M1-1	25ml
	PM50124M1-2	100ml
	PM50124M1-3	500ml
	PM50124M1-4	1L
	PM50124M1-5	5L
	PM50124M1-6	10L

Production date: See label Expiry date: 5 years, under proper storage conditions

#### Manufacturer: Sunresin New Materials Co. Ltd.

Add:No. 135, Jinye Rd, Xi'an Hi-tech Industrial Development Zone, Shaanxi, 710076, China <u>www.seplite.com</u> <u>www.sunresin.com</u> E-mail: info.lifescience@sunresin.com

All information set forth herein is for informational purposes only. This information is general descriptive(introductory) information of SUNRESIN and its related products, technologies and services. Neither shall constitute the guarantee of SUNRESIN and its affiliates to products, technologies and services in specific fields and specific application conditions results, unless otherwise expressly noted. SUNRESIN and its affiliates assumes no obligation or liability for the information in this document. Customer is responsible for judging whether the information is appropriate for Customer's concrete demand and are obliged to understand whether the use of these products, technologies and services is permitted by the laws and regulations of their countries and relevant regions. Unless expressly stated, no freedom from infringement of use any patent or trademark or intellectual property rights owned by SUNRESIN or its affiliated companies under this document is to be inferred.

