

1. Description

Seplife® LXMS 15Q is a polymeric resin for ion exchange chromatography characterized by strong chemical stability and high rigidity for use in polishing steps with or without organic solvents and at high flow rates

- Highly uniform particle size for high resolution
- Strong anion exchange resin based on a hydrophilic coated styrene/divinylbenzene backbone functionalized with quaternary ammonium (Q) groups
- Suitable for ion exchange chromatography in the separation of proteins, peptides, oligonucleotides and other small and medium size biomolecules
- High stability to CIP, organic solvents and pH
- Regulatory Support File (RSF) is available for Seplife LXMS 15Q

Seplife® LXMS 15Q is polymeric resin for ion exchange chromatography based on styrene/divinylbenzene functionalized with quaternary ammonium (Q) with a highly uniform particle size distribution (15 micron)

2. Properties

| Product | Seplife® LXMS 15Q |
|------------------------------------|--|
| Appearance | White to light yellow spherical beads |
| Type | Strong base anion - Quaternary amine |
| Matrix | Polystyrene/divinylbenzene |
| Ligand | Quaternary Amino |
| Ion exchange capacity (mmol/ml) | 0.07-0.10 |
| Particle size range (µm) | 13-17 |
| Typical pore size (Å) | 500 |
| pH stability | 2-12 (operation), 1-14 (CIP) |
| Chemical stability | Stable in commonly used aqueous ion exchange buffers |
| Flow rate* (cm/h) | Max 600 |
| Dynamic binding capacity** (mg/ml) | ≥ 55 |
| Shipped as | Slurry in 20% ethanol solution |

*Testing conditions: chromatography column 16mm×200mm; column bed height 100mm; Packing pressure 2 MPa; Mobile phase 0.5mol/L NaCl.

**Testing conditions: Chromatography column 8mm×100mm; Column bed height 100mm; Packing pressure 0.5 MPa; Mobile phase: 20mM Tris-HCl, 50mM NaCl, pH 8.0; Sample BSA; Retention time 2 min.

3. Instructions

3.1 Column loading

The best column packing results are obtained with the resin slurry in 0.5M NaCl at a concentration of 60%-70%.

- Mix the resins to form a homogenate and measure the desired mass or volume, about 1.2 times the column volume.
- Replace 20% ethanol with 0.5 M NaCl solution and equilibrate overnight.
- Before loading the column, adjust the slurry concentration to 65-70% with 0.5 M NaCl solution; pour the resin slurry into the chromatography column in one go, allow to sediment by gravity and mark the height after sedimentation.
- Load the distributor and adjust the height so that the compression coefficient is 1.05-1.10; then pump 0.5M NaCl through the column and stabilize the bed with 1.5 to 2 times the highest working flow rate. (2 - 5 column volumes).
- Efficiency and asymmetry determinations can be performed as described below.

3.2 Column Efficiency Evaluation

The test method for column efficiency of the chromatography columns is as follows:

Mobile phase: 0.5 M NaCl solution

Linear flow rate: 100 cm/h

Sample: 2 M NaCl solution

Loading volume: 0.5-1 % of column volume

Detection: Conductivity detector

Acceptable column asymmetry factor is in the range 0.8-1.8.

3.3 Rinsing

The loaded columns should be rinsed with at least 5 BV of deionized water.

3.4 Equilibration

Equilibrate the column with an appropriate 5-10 column volume buffer until the conductivity and pH of the effluent remain unchanged.

3.5 Sample feeding

The solid sample can be prepared by dissolving in the equilibrium solution. Low-concentration sample solutions can be concentrated in advance as much as possible while too high concentration sample solution can be diluted with the equilibration solution. To avoid clogging of the column, samples should be processed by centrifugation or membrane filtration. The feed amount is calculated according to the capacity of the resin and the content of the target protein in the feed solution. Before loading, ensure the sample buffer is as consistent as possible with the equilibration solution.

3.6 Elution

After loading the sample, continue rinsing with equilibration buffer until the baseline is stable. Typically, the molecules retained on the chromatographic media are eluted by increasing the salt concentration or changing the pH to ensure repulsion of the charged species. Occasionally, a solution containing a small percentage organic solvent such as acetonitrile or ethanol is required to disrupt any hydrophobic interactions.

3.7 Regeneration and CIP

Regular CIP can prevent the enrichment of protein precipitation impurities in the column bed, and help to maintain the capacity and separation effect of the chromatography media. In general, specific CIP methods and the frequency of CIP need to be designed for each process according to the type of contaminants. The recommended regeneration and CIP method is as follows: Rinse with 5 BV 1-2 M NaCl followed by 5 BV of 0.5-1 M NaOH.

4. Storage

Sealed and stored at 4-30°C (storage solution is 20% ethanol) in a ventilated, dry and clean place, do not freeze.

5. Transportation

Avoid sunlight, rain, and heavy pressure during transportation. It is strictly forbidden to transport with toxic and hazardous materials.

6. Ordering information

| Product Name | References | Pack Size |
|----------------------|-----------------|-----------|
| Seplife® LXMS 15Q | PS40323X(15)1-1 | 25ml |
| | PS40323X(15)1-2 | 100ml |
| | PS40323X(15)1-3 | 500ml |
| | PS40323X(15)1-4 | 1L |
| | PS40323X(15)1-5 | 5L |
| | PS40323X(15)1-6 | 10L |

Production date: See label

Expiry date: 5 years, under proper storage conditions

Manufacturer: Sunresin New Materials Co. Ltd.

Add:No. 135, Jinye Rd, Xi'an Hi-tech Industrial Development Zone, Shaanxi, 710076, China

www.seplife.com www.sunresin.com

E-mail: info.lifescience@sunresin.com

All information set forth herein is for informational purposes only. This information is general descriptive(introductory) information of SUNRESIN and its related products, technologies and services. Neither shall constitute the guarantee of SUNRESIN and its affiliates to products, technologies and services in specific fields and specific application conditions results, unless otherwise expressly noted. SUNRESIN and its affiliates assumes no obligation or liability for the information in this document. Customer is responsible for judging whether the information is appropriate for Customer's concrete demand and are obliged to understand whether the use of these products, technologies and services is permitted by the laws and regulations of their countries and relevant regions. Unless expressly stated, no freedom from infringement of use any patent or trademark or intellectual property rights owned by SUNRESIN or its affiliated companies under this document is to be inferred.