
# Seplife<sup>®</sup> Unylinker-OH Resin



### 1. Description

Seplife<sup>®</sup> Unylinker-OH Resin is a macroporous polystyrene resin functionalized with Unylinker groups with uniform particle size and active groups distribution. The Seplife<sup>®</sup> Unylinker-OH resin offers 99% coupling efficiency, it is fully compatible with multiple chemistries as well as standard solid phase oligonucleotide synthesis (SPOS) conditions and standard DMT deprotection conditions.



### 2. Properties

Product	Seplife <sup>®</sup> Unylinker-OH resin		
Appearance	Off white or light yellow beads		
Туре	Macroporous resin functionalized with Unylinker		
	groups		
Ligand	Unylinker		
Matrix	Polystyrene		
Substitution (micromol/g)	125-175		
	175-225		
	225-275		
	275-325		
	325-375		
	Customized		
Average Pore Size (nm)	30-60		
Swelling in Toluene (ml/g)	4.0-7.0		
Dry volume (ml/g)	2.0-3.0		
Particle size range (micron)	50 - 100		
Shipped as	Dry		



## Seplife<sup>®</sup> Unylinker-OH Resin



### 3.Storage

Seplife<sup>®</sup> Unylinker-OH Resin that is not for immediate use should be stored in dry form in a closed container at 2-8 °C.

### 4. Transportation

Avoid sunlight, rain, and heavy pressure during transportation. It is strictly forbidden to transport together with toxic and hazardous materials. Product can be shipped at ambient temperature.

### 5.Ordering information

Product Name	References	Substitution (µmol/g)
Seplife® Unylinker-OH Resin	LXSSN1-045-01-01	125-175
	LXSSN1-045-01-02	175-225
	LXSSN1-045-01-03	225-275
	LXSSN1-045-01-04	275-325
	LXSSN1-045-01-05	325-375
	LXSSN1-045-01-06	Customized

Production date: See label Retest date: 1 year, under proper storage conditions

Manufacturer: Sunresin New Materials Co. Ltd. Add:No. 135, Jinye Rd, Xi'an Hitech Industrial Development Zone, Shaanxi, 710076, China <u>www.seplite.com</u> <u>www.sunresin.com</u>

E-mail: info.lifescience@sunresin.com

All information set forth herein is for informational purposes only. This information is general descriptive(introductory) information of SUNRESIN and its related products, technologies and services. Neither shall constitute the guarantee of SUNRESIN and its affiliates to products, technologies and services in specific fields and specific application conditions results, unless otherwise expressly noted. SUNRESIN and its affiliates assumes no obligation or liability for the information in this document. Customer is responsible for judging whether the information is appropriate for Customer's concrete demand and are obliged to understand whether the use of these products, technologies and services is permitted by the laws and regulations of their countries and relevant regions. Unless expressly stated, no freedom from infringement of use any patent or trademark or intellectual property rights owned by SUNRESIN or its affiliated companies under this document is to be inferred.

