Seplife® EMC7120S



Descriptions

Seplife[®] EMC7120S is a highly cross-linked amino acrylate carrier, functionalized with a long spacer (C6) and that has the appearance of white spherical beads. The particle size distribution is 100-300 micron. The resin is highly hydrophilic and has a porous structure combined with a very high surface area which is ideal for immobilization of enzymes. Amino resins require pre-activation with glutaraldehyde and then aldehyde groups form multipoint covalent binding between the enzyme and resin. All amino acrylate resins are designed to form very stable covalent linkages with different protein groups (amino, thiol, phenolic) via Schiff base in very mild experimental conditions of pH and temperature. Amino acrylates can be used in either stirred tank or bed reactor.

The high cross-linking of the resin ensures extremely high mechanical stability thus allowing the possibility to reuse the immobilized enzyme for many cycles.

Matrix	Polyacrylate	
Appearance	White opaque spherical beads	
Functional group	Amino (C6)	
Immobilized method	Covalent bond	
Moisture content (%)	50-60	
Shipping weight (g/ml)	0.68-0.78	
Particle size (µm)	100-300	
Pore size (Å)	200-400	
Surface area (m ² /g)	≧100	
Chemical stability	The resin is insoluble in water, acid, alkali and methanol, ethanol, acetone, toluene, n-heptane, DMSO and other organic solvents.	
Immobilization pH	4-10	

Physical and Chemical Characteristics

•Key features and Benefits

- •Easy immobilization through pre-activation with glutaraldehyde
- ·Good physical and chemical stability
- Long lifetime

Recommended Storage Conditions

Recommended temperature of storage (°C)	2-20
Shelf life	3 years

Precautions

Resins should be stored in sealed containers or bags where temperature was 2-20°C in dry conditions without exposure to direct sunlight.

Do not mix ion exchange or adsorbent resin with strong oxidizing agents; otherwise it will cause violent reactions.



-ADD: No. 135, Jinye Rd, Xi'an Hi-tech Industrial Development Zone, Shaanxi, 710076, China -Email: info.lifescience@sunresin.com website: <u>www.sunresin.com</u>



Product Data Sheet

Seplife® EMC7120S



In case of eyes contact with resins, rinse eyes immediately with plenty of water, and consult a specialist. Material and samples must be disposed according to local regulations. Dry polymers will expand when become wetted and may cause an exothermic reaction. Spilled materials on the floor can cause slippery conditions.

•Ordering Information

Product Name	Reference number	Packing size
	PM033S01	500g
	PM033S02	1kg
Seplife [®] EMC7120S	PM033S03	5kg
	PM033S04	10kg
	PM033S05	25kg

·SUNRESIN® and SEPLIFE® are registered trademarks of Sunresin New Materials Co. Ltd.

•For more information about SEPLIFE® resins, please contact SUNRESIN® directly.





